



User Manual for Hydrocolloid Dressing



Instruction Book

【Ingredients】 medical hot melt glue, sodium carboxymethyl cellulose (CMC), polyisobutene

【Intended Use】 Apply to the injury to cover, keep the wound moist, and absorb the exudates, contact with skin epidermis.

【Instructions】

1. Check the package before use. Wash the wound with saline solution, and dry the skin around the skin wound gently .
2. Select the suitable size of the product (When necessary to cut the shape of the products), place the product on the wound completely.
3. Use your fingers to lightly press the edges a few times to reinforce the fixation.

【Contraindications】

1. People who are allergic to Hydrocolloid materials are forbidden, cautiously use of allergic constitution.
2. It is not recommended for dry wound.
3. Be careful with the infection of the wound.

【Cautions】

1. This product is a disposable sterile product. It is strictly prohibited to reuse the product. Packaging damage is disabled.
2. Please use it within the validity period. Read the instructions carefully before use.
3. Before use of this product for refractory wound, should make a scientific analysis of wound types and causes and formulate corresponding treatment plans. Only this product can not be used as a refractory wound treatment scheme.
4. Replace with a new one when it turns white and has absorbed tissue fluid

【Storage conditions】 The products shall be stored in non-corrosive gas and well

ventilated environment with room temperature.



3 years



see package



see package



Jiangsu Ostup Medical Products Co., Ltd.

No.9, Caofang Road,Luhe Zone, Nanjing,211500,Jiangsu, P.R. China

Tel: 0086-25-57678895 Website: www.ostupmedical.com

EU

REP

Prolinx GmbH

Brehmstr. 56, 40239, Duesseldorf, Germany

Tel: 0049 211 3105 4698



: Manufacturer



: Authorized representative in European



: Date of Manufacture



: Lot No.



: Expiry Date



Sterilized using irradiation



Do not use if package is damaged



Storage temperature



Disposable use



0197