



User Manual for Silicone Foam Dressing



Instruction Book

【Components】 Silicone Foam dressings are composed of silicone adhesive layer, hydrophilic polyurethane foam absorption layer, polyurethane film backing layer and the release film.

【Intended Use】 Be used to the wound care, cover the wound surface, absorb and control the exudate of the wound surface through the foam absorption layer, which can absorb a small to a medium amount of exudate, and also can be used to prevent pressure injury.

【Directions for use】

1. Clean or disinfect the wound before usage.
2. Select the appropriate size of the dressing, which the dressing should be 2cm larger than the wound;
3. Check and open the package, place the dressing on the skin, and then gently press the whole dressing, smooth it from the center to the periphery.
4. Keep observing the absorption of the dressing, replace the dressing in time when the dressing absorb adequate of the exudate.
5. Lift up a corner of the dressing, remove it.

【Contraindication】 Use with caution if allergic reaction occurs.

【Storage condition】 Store in a cool, dry, ventilated place away from direct sunlight.

【Cautions and warnings】

1. The product is sterilized, do not use if the package is damaged.
2. Disposable product, do not reuse.

3. Stop using the product beyond the expiration date.
4. Suggesting usage time should not exceed 7 days. If the liquid permeation volume is large, it should be replaced in time according to the actual situation.
5. The product should be used under the direction of the doctor or the nurse.
6. If you experience skin irritation and allergy, please stop using and consult medical professionals immediately.
7. Do not use the product with an oxidizing agent such as a hypochlorous acid solution or a catalase.
8. Please follow the hospital waste disposal requirements and rules to discard the product. Pay attention to protecting the environment.



3 years



see package



see package

Jiangsu Ostup Medical Products Co., Ltd.

No.9, Caofang Road,Luhe Zone, Nanjing,211500,Jiangsu, P.R. China

Tel: 0086-25-57678895 Website: www.ostupmedical.com

: Manufacturer



: See the instructions before usage



: Date of Manufacture



: Lot No.



: Expiry Date



: Do not re-use



: sterilized using ethylene oxide



: Do not use if the package is damaged



: Keep dry



: Keep away from sunlight